

Drug-induced liver injury: prediction using cellular imaging

Purpose

To clarify if long time (14 days) treatment of drugs in a functionally stable human hepatocyte co-culture on Cell-able[®] increase the predictive performance for clinical drug-induced liver injury (DILI) over that of the short time (24 hours) treatment ^a.

^a Xu, et al., 2008



Protocol

- 3T3 Swiss Albino cells (8x10³ cells/well) were seeded as feeder cells on Cell-able[®] (96well plates).
- 2. The next day, human hepatocytes $(2x10^4 \text{ cells/well})$ were seeded.
- 3. Three days later, vehicle or test drugs (1-60*Cmax) were treated.
- 4. The hepatocytes were maintained for 14 days. The medium was replaced by the new medium containing vehicle or test drugs three times a week.
- 5. The cells were stained with DRAQ5[®] for nuclei, monochlorobimane (mBCl) for GSH, tetramethylrhodamine methyl ester (TMRM) for mitochondrial membrane potential (MMP).
- 6. The cells were taken images by ImageXpress MICRO.



Acetaminophen-induced hepatocyte injury cultured on Cell-able®

	Нер	Hepatocytes + feeder cells				Feeder cells only			
	PC	Nuclei	GSH	MMP	PC	Nuclei	GSH	MMP	
		DRAQ5		TMRM	Light	DRAQ5	mBCI	TMRM	
Vehicle					e = = = + + + + + + + + + + + + + + + +				
(0.1% DMSO)					· 原生的中,也。				
				1	4 9 4 4 4 4 G				
Acetaminophen					· · · · · · · · · · · · · · · · · · ·				
Cmay	1 1 5 8 9 9 9 9 9 9 9 9 9 9 9 9 9 9 9 9 9 9 9			p.o. 0.1					
Cmax					安美海路部分 化				
					10 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0				
10*Cmax	0 * * * * * * * * * * * * * * * *			00000000000000000000000000000000000000	5 年間は回転 5 年間は				
				3 8 4 3 3 8 9 9	1 2 4 2 4 1 • • 2 • 2 • 2				
					3 3 2 4 4 3				
30*Cmax				100000 III	·····································				
					1949 - 1997 1949 - 1969 1949 - 1969				
	The for								
60*Cmax					1 2 1 d 2				
				******	1 4 3 to (2)				
				2000000	· · · · · · · · ·				

The cells were exposed to acetaminophen for 14 days.

Summary of sensitivity and specificity in terms of prediction of clinical DILI

- Human hepatocyte co-culture on Cell-able[®] maintained the liver functions, such as CYPs activities and albumin secretion for 54 days.
- The imaged GSH and MMP were sensitive parameters for DILI.
- The functionally stable hepatocyte co-culture on Cell-able[®] enabled the long time (14 days) exposure of drugs to hepatocytes, which conferred high predictive performance for clinical DILI.
- In this study, the sensitivity was about two-times higher than that in the previous report (Xu, 2008) in terms of prediction of clinical DILI.

	TRP, 2013	Xu, 2008
Clinical DILI		
Sensitivity	77% (10/13)	38% (5/13)
Specificity	100% (5/5)	100% (5/5)
Clinical DILI + LTKB		
Sensitivity	90% (9/10)	50% (5/10)
Specificity	88% (7/8)	100% (8/8)

Sensitivity: predicted positive/positive in clinical DILI, specificity: predicted negative/negative in clinical DILI