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

薬物性肝毒性 (Drug induced liver injury, DILI)
の予測

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
Molecular
Devices

ImageXpress Micro

High content screening assay system

TOYO GOSEI

Acetaminophen

8

The cells were exposed to acetaminophen for 14 days.

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




- Drugs of no concern for DILI ; N
- Drugs of less concern for DILI ; +/-
- Drugs of most concern for DILI; P

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10

| ID | 判定 LTKB+ Clinical DILI | Xu, 2008, 1-d, Matrigel, 100x Cmax | Khetani, 2013, 16-d, co-culture, 1-100x Cmax | Human hepatocyte, 14-d, co-culture Cell-able®, 1-30x Cmax | | PXB-cells™, 14-d, co-culture Cell-able®, 1-60x Cmax | |
|-------------------|------------------------------|--|--|---|-------------|---|----------|
| Aspirin | Negative | Negative | Negative | Negative | | Negative | |
| Lidocaine | Negative | Negative | Positive | Negative | | Negative | |
| Prednisone | Negative | Negative | Negative | Negative | | Negative | |
| Bupropion HCl | +/- | Negative | NA | Negative | | Negative | |
| Cyclophosphamide | +/- | Negative | Positive | Positive | GSH, MMP | Positive | GSH, MMP |
| Imipramine HCl | +/- | Negative | Positive | Negative | | Negative | |
| Propranolol HCl | +/- | Negative | Negative | Negative | | Negative | |
| Nifedipine | +/- | Negative | Negative | Negative | | Negative | |
| Acetaminophen | Positive | Positive | Positive | Positive | GSH, MMP, N | Positive | GSH, MMP |
| Benzbromarone | Positive | Positive | Positive | Positive | GSH, MMP, N | Positive | GSH, MMP |
| Ciprofloxacin HCl | Positive | Negative | Positive | Positive | GSH, MMP | Positive | GSH, MMP |
| Clomiplamine HCl | Positive | Negative | Positive | Positive | GSH, MMP | Negative | |
| Diclofenac Na | Positive | Positive | Positive | Positive | GSH, MMP, N | Positive | GSH, MMP |
| Estrone | Positive | Negative | Negative | Negative | | Negative | |
| Isoniazid | Positive | Negative | Positive | Positive | GSH, MMP | Positive | GSH |
| Phenacetin | Positive | Positive | Positive | Positive | GSH, MMP, N | Positive | GSH, MMP |
| Troglitazone | Positive | Positive | Positive | Positive | GSH, MMP, N | Positive | GSH, MMP |
| Ketoconazol | Positive | NA | Positive | Positive | GSH, MMP, N | Positive | GSH, MMP |

1) Xu et al. (2008) Cellular Imaging Predictions of Clinical Drug-Induced Liver Injury. *Toxicological Sciences* 105, 97-105

2) Khetani et al. (2013) The Use of Micropatterned Co-cultures to Detect Compounds that Cause Drug induced Liver Injury in Humans. *Toxicological Sciences* 132, 107-117

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| | Cell-able [®] | | Xu, 2008 ¹⁾ | Khetani, 2013 ²⁾ |
|------|------------------------|------------------------|------------------------|-----------------------------|
| | Human hepatocyte | PXB-cells [™] | Human hepatocyte | Human hepatocyte |
| 特異性* | 100% | 100% | 100% | 67% |
| 感度** | 90% | 80% | 56% | 90% |

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